




Ti·MESH CLINICAL CASES


수술 전 골결손부 방사선 사진



Ti·MESH를 적용한 골증대술



수술 5개월 후 치조골 방사선 사진




5개월 동안 증대된 치조골

PRODUCT FEATURES

01. 최소한의 고정으로 간편하게 사용
02. Window가 형성되어 Ti·MESH를 고정 후 골이식이 가능함
03. 타사제품과 비교하여 Pore이 작아 Ti·MESH 제거하기가 수월함
04. Ti·MESH 제거 후에도 판막상태가 양호하여 추가 골증대술 가능

MULTI UNIT




Product	Ref.	L(mm)	W(mm)	D(φ)	T(mm)
D Type-S	TM362210	36.0	22.0	1.7	0.2
D Type-M	TM363310	36.0	33.0	1.7	0.2

Bone Augmentation



SINGLE UNIT




Product	Ref.	L(mm)	W(mm)	D(φ)	T(mm)
A Type-S	TM36110	36.0	10.0	1.7	0.2
A Type-M	TM3612	36.0	12.0	1.7	0.2

Vertical Bone Augmentation





Horizontal Bone Augmentation



Ti·MESH CLINICAL CASES


Preoperative radiographs of bone defects






Radiographs of Expanded alveolar bone in 5 months.

Osteopathic surgery by using Ti·MESH






Expanded alveolar bone in 5 months

PRODUCT FEATURES

- 01. Utilized as minimum fixing
- 02. Osseintegration after fixing
Ti·MESH is available due to window in mesh.
- 03. Simple removing Ti·MESH due to small pore
- 04. Additional Osteopathic surgery is available due to good condition after eliminating Ti·MESH

MULTI UNIT



Product	Ref.	L(mm)	W(mm)	D(ø)	T(mm)
D Type-S	TM362210	36.0	22.0	1.7	0.2
D Type-M	TM363310		33.0		


Bone Augmentation







SINGLE UNIT



Product	Ref.	L(mm)	W(mm)	D(ø)	T(mm)
A Type-S	TM36110	36.0	10.0	1.7	0.2
A Type-M	TM3612		12.0		

Vertical Bone Augmentation







Horizontal Bone Augmentation





